

EGR3 Antibody, HRP conjugated



PACO58625

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable transcription factor involved in muscle spindle development.

Gene ID:

EGR3

Uniprot

Q06889

Synonyms:

Early growth response protein 3 (EGR-3) (Zinc finger protein pilot), EGR3, PILOT

Immunogen:

Recombinant Human Early growth response protein 3 protein (201-281AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A