

LECT2 Antibody, Biotin conjugated



PACO58812

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has a neutrophil chemotactic activity. Also a positive regulator of chondrocyte proliferation. Does not show metalloendopeptidase activity.

Gene ID:

LECT2

Uniprot

O14960

Synonyms:

Leukocyte cell-derived chemotaxin-2 (LECT-2) (hLECT2), LECT2

Immunogen:

Recombinant Human Leukocyte cell-derived chemotaxin-2 protein (24-151AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A