

RANBP10 Antibody, HRP conjugated



PACO59166

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a guanine nucleotide exchange factor (GEF) for RAN GTPase. May play an essential role in hemostasis and in maintaining microtubule dynamics with respect to both platelet shape and function. May act as an adapter protein to couple membrane receptors to intracellular signaling pathways. Enhances dihydrotestosterone-induced transactivation activity of AR, as well as dexamethasone-induced transactivation activity of NR3C1, but does not affect estrogen-induced transactivation. In contrast to RANBP9, does not interact with Sos and does not activate the Ras pathway.

Gene ID:

RANBP10

Uniprot

Q6VN20

Synonyms:

Ran-binding protein 10 (RanBP10), RANBP10, KIAA1464

Immunogen:

Recombinant Human Ran-binding protein 10 protein (3-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A