

VRK2 Antibody, Biotin conjugated



PACO59176

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

VRK2

Uniprot

Q86Y07

Synonyms:

Serine/threonine-protein kinase VRK2 (EC 2.7.11.1) (Vaccinia-related kinase 2), VRK2

Immunogen:

Recombinant Human Serine/threonine-protein kinase VRK2 protein (317-472AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A