

RASSF3 Antibody, Biotin conjugated



PACO59204

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, plasma membrane, identical protein binding.

Gene ID:

RASSF3

Uniprot

Q86WH2

Synonyms:

Ras association domain-containing protein 3, RASSF3

Immunogen:

Recombinant Human Ras association domain-containing protein 3 protein (2-109AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A