

FAM71D Antibody




PACO59233

Product Information

Size:	Protein Background:
50ug	nucleus.
Reactivity:	Gene ID:
Human	FAM71D
Source:	Uniprot
Rabbit	Q8N9W8
Isotype:	Synonyms:
IgG	Protein FAM71D, FAM71D, C14orf54
Applications:	Immunogen:
ELISA, WB	Recombinant Human Protein FAM71D protein (220-398AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:500-1:5000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: MCF-7 whole cell lysate. All lanes: FAM71D antibody at 4.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 48, 45 kDa. Observed band size: 48 kDa.