

LYRM4 Antibody, HRP conjugated



PACO59426

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for nuclear and mitochondrial iron-sulfur protein biosynthesis.

Gene ID:

LYRM4

Uniprot

Q9HD34

Synonyms:

LYR motif-containing protein 4, LYRM4, C6orf149 ISD11

Immunogen:

Recombinant Human LYR motif-containing protein 4 protein (1-91AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A