

HES7 Antibody, FITC conjugated



PACO59515

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional repressor. Represses transcription from both N box- and E box-containing promoters. May with HES1, cooperatively regulate somite formation in the presomitic mesoderm (PSM). May function as a segmentation clock, which is essential for coordinated somite segmentation.

Gene ID:

HES7

Uniprot

Q9BYE0

Synonyms:

Transcription factor HES-7 (hHes7) (Class B basic helix-loop-helix protein 37) (bHLHb37) (Hair and enhancer of split 7) (bHLH factor Hes7), HES7, BHLHB37

Immunogen:

Recombinant Human Transcription factor HES-7 protein (127-192AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A