

SLC45A2 Antibody, Biotin conjugated



PACO59548

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Melanocyte differentiation antigen. May transport substances required for melanin biosynthesis.

Gene ID:

SLC45A2

Uniprot

Q9UMX9

Synonyms:

Membrane-associated transporter protein (Melanoma antigen AIM1) (Protein AIM-1) (Solute carrier family 45 member 2), SLC45A2, AIM1 MATP

Immunogen:

Recombinant Human Membrane-associated transporter protein (1-46AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A