

MKL2 Antibody, FITC conjugated



PACO59563

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a transcriptional coactivator of serum response factor (SRF). Required for skeletal myogenic differentiation.

Gene ID:

MKL2

Uniprot

Q9ULH7

Synonyms:

MKL/myocardin-like protein 2 (Megakaryoblastic leukemia 2) (Myocardin-related transcription factor B) (MRTF-B), MKL2, KIAA1243 MRTFB

Immunogen:

Recombinant Human MKL/myocardin-like protein 2 protein (164-313AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A