

PLAGL2 Antibody, Biotin conjugated



PACO59568

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Shows weak transcriptional activatory activity.

Gene ID:

PLAGL2

Uniprot

Q9UPG8

Synonyms:

Zinc finger protein PLAGL2 (Pleiomorphic adenoma-like protein 2), PLAGL2, KIAA0198

Immunogen:

Recombinant Human Zinc finger protein PLAGL2 protein (324-469AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A