

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a key role in retinoic acid, metabolism. Acts on retinoids, including all-trans-retinoic acid, (RA) and its stereoisomer 9-cis-RA. Capable of 4-hydroxylation. Responsible for generation of several hydroxylated forms of RA, including 4-OH-RA and 4-oxo-RA.

**Gene ID:**

cyp26a1

**Uniprot**

P79739

**Synonyms:**

Cytochrome P450 26A1 (EC 1.14.13) (Cytochrome P450RAI) (Retinoic acid, 4-hydroxylase) (Retinoic acid, metabolizing cytochrome), cyp26a1, cyp26

**Immunogen:**

Recombinant Zebrafish Cytochrome P450 26A1 protein (244-371AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A