

KRIT1 Antibody, HRP conjugated



PACO59769

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KRIT1

Uniprot

O00522

Synonyms:

Krev interaction trapped protein 1 (Krev interaction trapped 1) (Cerebral cavernous malformations 1 protein), KRIT1, CCM1

Immunogen:

Recombinant Human Krev interaction trapped protein 1 protein (471-585AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A