

PRB1 Antibody



PACO59864

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

PRB1

Uniprot

P04280

Synonyms:

Basic salivary proline-rich protein 1 (Salivary proline-rich protein) [Cleaved into: Proline-rich peptide II-2; Basic peptide IB-6; Peptide P-H], PRB1


Immunogen:

Recombinant Human Basic salivary proline-rich protein 1 protein (17-92AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: K562 whole cell lysate. All lanes: PRB1 antibody at 3.16 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.