

C14orf177 Antibody, HRP conjugated



PACO60065

Product Information

Size:

50ug

Protein Background:**Reactivity:**

Human

Gene ID:

C14orf177

Source:

Rabbit

Uniprot

Q52M58

Isotype:

IgG

Synonyms:

Putative uncharacterized protein C14orf177, C14orf177

Applications:

ELISA

Immunogen:

Recombinant Human Putative uncharacterized protein C14orf177 protein (1-125AA).

Recommended dilutions:**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A