

PDILT Antibody



PACO60220

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Probable redox-inactive chaperone involved in spermatogenesis.

Gene ID:

PDILT

Uniprot

Q8N807

Synonyms:

Protein disulfide-isomerase-like protein of the testis, PDILT


Immunogen:

Recombinant Human Protein disulfide-isomerase-like protein of the testis protein (494-584AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: A549 whole cell lysate, Colo320 whole cell lysate. All lanes: PDILT antibody at 4.63 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 67 kDa. Observed band size: 67 kDa.