

NKX2-3 Antibody



PACO60288

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Transcription factor.

Gene ID:

NKX2-3

Uniprot

Q8TAU0

Synonyms:

Homeobox protein Nkx-2.3 (Homeobox protein NK-2 homolog C), NKX2-3, NKX23, NKX2C

Immunogen:

Recombinant Human Homeobox protein Nkx-2.3 protein (52-148AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: HeLa whole cell lysate. All lanes: NKX2-3 antibody at 3.8 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 39 kDa. Observed band size: 39 kDa.