

TMEM87B Antibody, FITC conjugated



PACO60851

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in retrograde transport from endosomes to the trans-Golgi network (TGN).

Gene ID:

TMEM87B

Uniprot

Q96K49

Synonyms:

Transmembrane protein 87B, TMEM87B

Immunogen:

Recombinant Human Transmembrane protein 87B protein (451-554AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A