

FMN1 Antibody, HRP conjugated



PACO60986

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in the formation of adherens junction and the polymerization of linear actin cables.

Gene ID:

FMN1

Uniprot

Q68DA7

Synonyms:

Formin-1 (Limb deformity protein homolog), FMN1, FMN LD

Immunogen:

Recombinant Human Formin-1 protein (354-487AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A