

UNC5D Antibody, FITC conjugated



PACO61015

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for the netrin NTN4 that promotes neuronal cell survival. Plays a role in cell-cell adhesion and cell guidance. Receptor for netrin involved in cell migration. Plays a role in axon guidance by mediating axon repulsion of neuronal growth cones in the developing nervous system upon ligand binding. May play a role in apoptosis in response to DNA damage. It also acts as a dependence receptor required for apoptosis induction when not associated with netrin ligand. Mediates cell-cell adhesion via its interaction with FLRT3 on an adjacent cell.

Gene ID:

UNC5D

Uniprot

Q6UXZ4

Synonyms:

Netrin receptor UNC5D (Protein unc-5 homolog 4) (Protein unc-5 homolog D), UNC5D, KIAA1777 UNC5H4

Immunogen:

Recombinant Human Netrin receptor UNC5D protein (90-212AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A