
Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Plays a role in myogenic differentiation by regulating myog levels. May bind to the 3'-UTR of myog mRNA and regulate its stability.

Gene ID:

rbm24

Uniprot

Q76LC6

Synonyms:

RNA-binding protein 24 (RNA-binding motif protein 24) (RNA-binding protein SEB-4), rbm24


Immunogen:

Recombinant Zebrafish RNA-binding protein 24 protein (2-78AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: rbm24 antibody at 4.3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 41 kDa. Observed band size: 41 kDa.