

# IFT52 Antibody, HRP conjugated



PACO61266

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in ciliogenesis as part of a complex involved in intraflagellar transport (IFT), the bi-directional movement of particles required for the assembly, maintenance and functioning of primary cilia. Required for the anterograde transport of IFT88.

**Gene ID:**

IFT52

**Uniprot**

Q9Y366

**Synonyms:**

Intraflagellar transport protein 52 homolog (Protein NGD5 homolog), IFT52, C20orf9 NGD5

**Immunogen:**

Recombinant Human Intraflagellar transport protein 52 homolog protein (325-434AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A