

ZNF37A Antibody



PACO61314

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF37A

Uniprot

P17032

Synonyms:

Zinc finger protein 37A, Zinc finger protein KOX21, ZNF37A, KOX21, ZNF37


Immunogen:

Recombinant Human Zinc finger protein 37A protein (263-386AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: U87 whole cell lysate, K562 whole cell lysate, 293T whole cell lysate. All lanes: ZNF37A antibody at 2.6 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 66 kDa. Observed band size: 66 kDa.