

GALNT9 Antibody, HRP conjugated



PACO61443

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the initial reaction in O-linked oligosaccharide biosynthesis, the transfer of an N-acetyl-D-galactosamine residue to a serine or threonine residue on the protein receptor. Does not glycosylate apomucin or SDC3.

Gene ID:

GALNT9

Uniprot

Q9HCQ5

Synonyms:

Polypeptide N-acetylgalactosaminyltransferase 9, Polypeptide GalNAc transferase 9, GalNAc-T9, pp-GaNTase 9, Protein-UDP acetylgalactosaminyltransferase 9, polypeptide N-acetylgalactosaminyltransferase 9, GALNT9

Immunogen:

Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 9 protein (184-338AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A