

PACO61450

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May function as a transcription factor.

Gene ID:

mycn

Uniprot

Q9PSJ0

Synonyms:

N-myc protein, zN-Myc, mycn, nmyc1


Immunogen:

Recombinant Zebrafish N-myc protein (1-93AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: mycn antibody at 2.2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 41 kDa. Observed band size: 48 kDa.