

HES2 Antibody



PACO61494

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Transcriptional repressor of genes that require a bHLH protein for their transcription.

Gene ID:

HES2

Uniprot

Q9Y543

Synonyms:

Transcription factor HES-2, Class B basic helix-loop-helix protein 40, bHLHb40, Hairly and enhancer of split 2, HES2, BHLHB40


Immunogen:

Recombinant Human Transcription factor HES-2 protein (1-76AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of PC-3 cells with PACO61494 at 1:66, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).