

HRG Antibody, HRP conjugated



PACO61543

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	HRG
Source:	Uniprot
Rabbit	P04196
Isotype:	Synonyms:
IgG	Histidine-rich glycoprotein, Histidine-proline-rich glycoprotein, HPRG, HRG
Applications:	Immunogen:
ELISA	Recombinant Human Histidine-rich glycoprotein protein (412-511AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A