

TTC39B Antibody, FITC conjugated



PACO61716

Product Information

Size:	Protein Background:
50ul	
Reactivity:	Gene ID:
Human	TTC39B
Source:	Uniprot
Rabbit	Q5VTQ0
Isotype:	Synonyms:
IgG	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B), TTC39B, C9orf52
Applications:	Immunogen:
ELISA	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A