

P Antibody



PACO61862

Product Information

Size:

50ul

Reactivity:

Vesicular stomatitis Indiana virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May play a role in viral pathogenesis or transmission by insects vectors.

Gene ID:

P

Uniprot

Q86132

Synonyms:

Protein C', P


Immunogen:

Recombinant Vesicular stomatitis Indiana virus Protein C' protein (1-67AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: P antibody at 2.7 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 14 kDa. Observed band size: 14 kDa.