

vcp Antibody, FITC conjugated



PACO61872

Product Information

Size:

50ul

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

vcp

Uniprot

Q7ZU99

Synonyms:

Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (Protein CDC48) (Valosin-containing protein) (VCP), vcp, cdc48

Immunogen:

Recombinant Zebrafish Transitional endoplasmic reticulum ATPase protein (704-806AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A