

# vcp Antibody, Biotin conjugated



PACO61873

---

## Product Information

**Size:**

50ul

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

vcp

**Uniprot**

Q7ZU99

**Synonyms:**

Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (Protein CDC48) (Valosin-containing protein) (VCP), vcp, cdc48

**Immunogen:**

Recombinant Zebrafish Transitional endoplasmic reticulum ATPase protein (704-806AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A