

CRH Antibody, Biotin conjugated



PACO61922

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Hormone regulating the release of corticotropin from pituitary gland. Induces NLRP6 in intestinal epithelial cells, hence may influence gut microbiota profile.

Gene ID:

CRH

Uniprot

P06850

Synonyms:

Corticoliberin (Corticotropin-releasing factor) (CRF) (Corticotropin-releasing hormone), CRH

Immunogen:

Peptide sequence from Human Corticoliberin protein (29-45AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A