

TGFBR3L Antibody



PACO62019

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

TGFBR3L

Uniprot

H3BV60

Synonyms:

Transforming growth factor-beta receptor type 3-like protein (TGF-beta receptor type-3-like protein) (TGFR-3L) (Transforming growth factor-beta receptor type III-like protein) (TGF-beta receptor type III-like protein), TGFBR3L


Immunogen:

Peptide sequence from Human Transforming growth factor-β receptor type 3-like protein (242-258AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of PC3 cells with PACO62019 at 1:50, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).