

# TGFBR3L Antibody, FITC conjugated



PACO62021

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TGFBR3L

**Uniprot**

H3BV60

**Synonyms:**

Transforming growth factor-beta receptor type 3-like protein (TGF-beta receptor type-3-like protein) (TGFR-3L) (Transforming growth factor-beta receptor type III-like protein) (TGF-beta receptor type III-like protein), TGFBR3L

**Immunogen:**

Peptide sequence from Human Transforming growth factor-&beta; receptor type 3-like protein (242-258AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A