

CHN48 Antibody



PACO62031

Product Information

Size:

50ul

Reactivity:

Nicotiana tabacum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Defense against chitin-containing fungal pathogens.

Gene ID:

CHN48

Uniprot

P08252

Synonyms:

Endochitinase A (CHN-A) (EC 3.2.1.14), CHN48


Immunogen:

Recombinant Nicotiana tabacum Endochitinase A protein (24-322AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Cucurbita moschata leaf tissue, Solanum tuberosum leaf tissue, Solanum lycopersicum leaf tissue. All lanes: CHN48 antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 36 kDa. Observed band size: 36 kDa.