

CAT1 Antibody, FITC conjugated



PACO62073

Product Information

Size:

50ul

Reactivity:

Arabidopsis thaliana

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Occurs in almost all aerobically respiring organisms and serves to protect cells from the toxic effects of hydrogen peroxide.

Gene ID:

CAT1

Uniprot

Q96528

Synonyms:

Catalase-1 (EC 1.11.1.6), CAT1

Immunogen:

Recombinant Arabidopsis thaliana Catalase-1 protein (1-492AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A