

ITGB1BP1 Antibody, Biotin conjugated



PACO62318

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

ITGB1BP1

Uniprot

O14713

Synonyms:

Integrin beta-1-binding protein 1 (Integrin cytoplasmic domain-associated protein 1) (ICAP-1), ITGB1BP1, ICAP1

Immunogen:

Recombinant Human Integrin beta-1-binding protein 1 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A