

# PCSK9 Antibody, FITC conjugated



PACO62333

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCSK9

**Uniprot**

Q8NBP7

**Synonyms:**

Proprotein convertase subtilisin/kexin type 9 (EC 3.4.21) (Neural apoptosis-regulated convertase 1) (NARC-1) (Proprotein convertase 9) (PC9) (Subtilisin/kexin-like protease PC9), PCSK9, NARC1

**Immunogen:**

Recombinant Human Proprotein convertase subtilisin/kexin type 9 protein (153-380AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A