

TBKBP1 Antibody



PACO62387

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Adapter protein which constitutively binds TBK1 and IKKε playing a role in antiviral innate immunity.

Gene ID:

TBKBP1

Uniprot

A7MCY6

Synonyms:

TANK-binding kinase 1-binding protein 1 (TBK1-binding protein 1), TBKBP1, KIAA0775
SINTBAD


Immunogen:

Recombinant Human TANK-binding kinase 1-binding protein 1 protein (327-427AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of HeLa cells with PACO62387 at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).