

BCL2L10 Antibody, Biotin conjugated



PACO62626

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes cell survival. Suppresses apoptosis induced by BAX but not BAK.

Gene ID:

BCL2L10

Uniprot

Q9HD36

Synonyms:

Bcl-2-like protein 10 (Bcl2-L-10) (Anti-apoptotic protein NrH) (Apoptosis regulator Bcl-B), BCL2L10, BCLB

Immunogen:

Recombinant Human Bcl-2-like protein 10 protein (19-127AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A