

TBC1D25 Antibody, FITC conjugated



PACO62761

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a GTPase-activating protein specific for RAB33B. Involved in the regulation of autophagosome maturation, the process in which autophagosomes fuse with endosomes and lysosomes.

Gene ID:

TBC1D25

Uniprot

Q3MII6

Synonyms:

TBC1 domain family member 25, TBC1D25, OATL1

Immunogen:

Recombinant Human TBC1 domain family member 25 protein (1-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A