

PIP4K2C Antibody, Biotin conjugated



PACO62770

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play an important role in the production of Phosphatidylinositol bisphosphate (PIP2), in the endoplasmic reticulum.

Gene ID:

PIP4K2C

Uniprot

Q8TBX8

Synonyms:

Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma (EC 2.7.1.149)
(Phosphatidylinositol 5-phosphate 4-kinase type II gamma) (PI(5)P 4-kinase type II gamma) (PIP4KII-gamma), PIP4K2C, PIP5K2C

Immunogen:

Recombinant Human Phosphatidylinositol 5-phosphate 4-kinase type-2 γ protein (127-408AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A