

PXN Antibody, Biotin conjugated



PACO62966

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion).

Gene ID:

PXN

Uniprot

P49023

Synonyms:

Paxillin, PXN

Immunogen:

Peptide sequence from Human Paxillin protein (28-45AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A