

RASSF2 Antibody, FITC conjugated



PACO62973

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Potential tumor suppressor. Acts as a KRAS-specific effector protein. May promote apoptosis and cell cycle arrest. Stabilizes STK3/MST2 by protecting it from proteasomal degradation.

Gene ID:

RASSF2

Uniprot

P50749

Synonyms:

Ras association domain-containing protein 2, RASSF2, CENP-34 KIAA0168

Immunogen:

Recombinant Human Ras association domain-containing protein 2 protein (69-326AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A