

TFPI2 Antibody, HRP conjugated



PACO63000

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the regulation of plasmin-mediated matrix remodeling. Inhibits trypsin, plasmin, factor VIIa/tissue factor and weakly factor Xa. Has no effect on thrombin.

Gene ID:

TFPI2

Uniprot

P48307

Synonyms:

Tissue factor pathway inhibitor 2 (TFPI-2) (Placental protein 5) (PP5), TFPI2

Immunogen:

Recombinant Human Tissue factor pathway inhibitor 2 protein (23-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A