

# MDM2 Antibody, HRP conjugated



PACO63016

---

## Product Information

|                               |   |
|-------------------------------|---|
| <b>Size:</b>                  | <b>Protein Background:</b>  |
| 50ul                          |   |
| <b>Reactivity:</b>            | <b>Gene ID:</b>   |
| Human                         | MDM2  |
| <b>Source:</b>                | <b>Uniprot</b>  |
| Rabbit                        | Q00987  |
| <b>Isotype:</b>               | <b>Synonyms:</b>  |
| IgG                           | E3 ubiquitin-protein ligase Mdm2 (EC 2.3.2.27) (Double minute 2 protein) (Hdm2) (Oncoprotein Mdm2) (RING-type E3 ubiquitin transferase Mdm2) (p53-binding protein Mdm2), MDM2 |
| <b>Applications:</b>          | <b>Immunogen:</b>   |
| ELISA                         | Peptide sequence from Human E3 ubiquitin-protein ligase Mdm2 protein (374-391AA).   |
| <b>Recommended dilutions:</b> | <b>Storage:</b>   |
|                               | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4  |

## Product Images

---

N/A

N/A