

GPC3 Antibody, Biotin conjugated



PACO63022

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cell surface proteoglycan that bears heparan sulfate. Inhibits the dipeptidyl peptidase activity of DPP4. May be involved in the suppression/modulation of growth in the predominantly mesodermal tissues and organs. May play a role in the modulation of IGF2 interactions with its receptor and thereby modulate its function. May regulate growth and tumor predisposition.

Gene ID:

GPC3

Uniprot

P51654

Synonyms:

Glypican-3 (GTR2-2) (Intestinal protein OCI-5) (MXR7) [Cleaved into: Secreted glypican-3], GPC3, OCI5

Immunogen:

Recombinant Human Glypican-3 protein (346-550AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A