

fkpA Antibody, Biotin conjugated



PACO63062

Product Information

Size:

50ul

Reactivity:

Aeromonas hydrophila

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. FkpA probably acts in the folding of extracytoplasmic proteins.

Gene ID:

fkpA

Uniprot

O08437

Synonyms:

FKBP-type peptidyl-prolyl cis-trans isomerase FkpA (PPIase) (EC 5.2.1.8) (Rotamase), fkpA

Immunogen:

Recombinant Aeromonas hydrophila FKBP-type peptidyl-prolyl cis-trans isomerase FkpA protein (21-268AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A