

NKAPL Antibody



PACO63099

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Transcriptional repressor of Notch-mediated signaling. Required for spermatogenesis.

Gene ID:

NKAPL

Uniprot

Q5M9Q1

Synonyms:

NKAP-like protein, NKAPL, C6orf194


Immunogen:

Peptide sequence from Human NKAP-like protein (114-123AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: A549 whole cell lysate, K562 whole cell lysate. All lanes: NKAPL antibody at 1:1000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 47 kDa. Observed band size: 47 kDa.