

ENHO Antibody, Biotin conjugated



PACO63126

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the regulation of glucose homeostasis and lipid metabolism.

Gene ID:

ENHO

Uniprot

Q6UWT2

Synonyms:

Adropin (Energy homeostasis-associated protein), ENHO, C9orf165

Immunogen:

Peptide sequence from Human Adropin protein (37-54AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A